

# Custom Ear Piece Order Form

Marvel™

## Customer Information

Ship To Account Number: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_  
 Bill To Account Number: \_\_\_\_\_  
 Third Party Patient Number: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 Purchase Order Number: \_\_\_\_\_  
 Contact Name: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Email: \_\_\_\_\_

## Patient Information

Last Name: \_\_\_\_\_  
 First Name: \_\_\_\_\_  
 Age: \_\_\_\_\_ Gender: \_\_\_\_\_

## Audiogram (Required for AOV):

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

## HI Warranty/Rush Options

2<sup>nd</sup> year  3<sup>rd</sup> year  4<sup>th</sup> year /  24-hour service (\$59)

## Technology Level

M90 (Premium)  M70 (Advanced)  M50 (Standard)  M30 (Essential)

## Hearing Device Selection – RIC

Qty	Instrument Type	Instrument Style			
		312	312T	13T	R
_____	Phonak Audéo™ M	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Specify color choice: \_\_\_\_\_

## Custom Ear Piece Selection – RIC

Side	Shell Styles	Power Level			
L	R	S	M (Std.)	P	UP
		46/111	50/114	58/122	67/130
<input type="checkbox"/>	<input type="checkbox"/> RIC SecureFit (Solid Silicone) [062-0007-01]	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A
<input type="checkbox"/>	<input type="checkbox"/> SlimTip (Hollow Acrylic) [062-0006-01]	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A
<input type="checkbox"/>	<input type="checkbox"/> SlimTip Titanium (Hollow) [062-0015-01]	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A
<input type="checkbox"/>	<input type="checkbox"/> cShell (Acrylic) [062-0008-01]	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Receiver Length (0-3): L \_\_\_\_\_ R \_\_\_\_\_

## Colors



## Ear Piece Options

### Shell Color Options

	Pink		Tan		Cocoa		Brown		Transp.		Titanium		Other <sup>1</sup>	
	L	R	L	R	L	R	L	R	L	R	L	R	L	R
RIC SecureFit (Silicone)	N/A	N/A	N/A	N/A	N/A	N/A	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A	N/A	N/A	N/A
SlimTip (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	<input type="checkbox"/>
SlimTip Titanium	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	<input type="checkbox"/>	<input type="checkbox"/>	N/A	N/A
cShell (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A	N/A	<input type="checkbox"/>	<input type="checkbox"/>

Specify color choice: \_\_\_\_\_

<sup>1</sup> Other color options: Red, Blue, White.

## Ear Piece Options, continued

### Faceplate Color Options

	Pink		Tan		Cocoa		Brown	
	L	R	L	R	L	R	L	R
cShell (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	(Std.)							

### Shell Options\*

	Canal Lock <sup>2</sup>		Skeleton Lock <sup>2</sup>		Helix Lock <sup>2</sup>		Extra Large	
	L	R	L	R	L	R	L	R
RIC SecureFit (Silicone)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A		N/A	
SlimTip (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A	
cShell (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

\*N/A for Titanium Shells

### Venting Options

	AOV (Audiogram required)		Pressure Vent		2.0mm SAV		2.5mm SAV		3.0mm SAV		Customized Large Vent		Cavity Vent		No Vent	
	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R
RIC SecureFit (Silicone)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A		<input type="checkbox"/>	<input type="checkbox"/>
SlimTip (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
SlimTip Titanium	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
cShell (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A		<input type="checkbox"/>	<input type="checkbox"/>
	(Std.)															

### Wax Options\*

	No Wax Prevention		Cerustop		Extended Receiver Tube		Wax Spring		Smart Guard <sup>3</sup> (Chargeable)	
	L	R	L	R	L	R	L	R	L	R
SlimTip (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A	
SlimTip Titanium	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	N/A	
cShell (Acrylic)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	(Std.)									

### Removal Line Color Options

	L	R
Transparent	<input type="checkbox"/>	<input type="checkbox"/>
Tan/Pink	<input type="checkbox"/>	<input type="checkbox"/>
Brown/Cocoa	<input type="checkbox"/>	<input type="checkbox"/>
No Removal Line	<input type="checkbox"/>	<input type="checkbox"/>
	(Std.)	

\*N/A for RIC SecureFit

<sup>2</sup> The shell and lock will be the same color.

<sup>3</sup> Compatible with S and M receivers only.

### Accessories

- TV Connector [076-3006-0611]
- Mini Charger [075-3009-11]
- Power Pack [071-0003]



### Preferences

- If necessary, may we change the following:  Please Call
- Yes      No
- UP to P OR P to M OR M to S if audiogram permits
- Build larger if required

### Measuring Tips

1. Make sure you are using the Phonak measuring tool
2. Make sure Receiver length measuring tool is parallel to the floor (straight on top of the ear)
3. Do not apply pressure to the tool (let it sit on the ear)
4. Use block colors to determine size (if on the line, use the smaller size)
5. Measure to the top of the ear canal opening



RIC SecureFit (Silicone)



SlimTip (Acrylic)



SlimTip Titanium



cShell (Acrylic)



cShell XL (Acrylic HS)

All of our products, including custom products and spare parts, can be ordered online in the Phonak store.

Sonova USA Inc. is not responsible and assumes no liability for any non-Phonak manufactured device or accessory sent by you to Phonak. Please ensure that you only include Phonak devices and accessories herein.