

Unitron Government Order Form

Group 3 (RIC)

unitron™

Step 1 Ship to

Ship to account number _____

Account name _____

Address _____

City _____ State _____ Zip code _____

Contact information

Date _____

Contact Name _____

Phone _____

Email address _____

Rush ☐ 24 Hour Rush ☐ 48 Hour Rush

Step 2 Patient information

First _____ Last 4 of SSN _____

Last _____

PO# _____

Patient audiogram

	250 Hz	500 Hz	1000 Hz	2000 Hz	3000 Hz	4000 Hz
Left						
Right						

Step 3 Hearing instrument selection

☐ L ☐ R Earmolds only (\$31.05)

North Platform \$375.71

Moxi™ Now Pro 050-5734	<input type="checkbox"/> L <input type="checkbox"/> R	Moxi Fit Pro 050-5694	<input type="checkbox"/> L <input type="checkbox"/> R
Moxi Kiss Pro 050-5684	<input type="checkbox"/> L <input type="checkbox"/> R	Moxi Dura Pro 050-5674	<input type="checkbox"/> L <input type="checkbox"/> R

Color (color samples on back)

- | | | |
|--|---|--|
| <input type="checkbox"/> Beige (01) | <input type="checkbox"/> Pewter (P7) | <input type="checkbox"/> Espresso boost (S3) |
| <input type="checkbox"/> Amber (P2) | <input type="checkbox"/> Charcoal (P8) | <input type="checkbox"/> Pewter shine (S5) |
| <input type="checkbox"/> Espresso (P4) | <input type="checkbox"/> Cinnamon (Q9) | <input type="checkbox"/> Sand storm (S6) |
| <input type="checkbox"/> Platinum (P6) | <input type="checkbox"/> Amber suede (S2) | <input type="checkbox"/> Teal blast (S7) |

Warranty (applicable on instruments and remotes only)

3 Year Warranty (standard, no charge)

Step 4 Wireless accessories

- | | |
|--|---|
| <input type="checkbox"/> Remote control 2 \$103.50
076-5031-0000070 | <input type="checkbox"/> uDirect™ 3 \$199.00 076-5033-T809070 |
| <input type="checkbox"/> uStream \$199.00 076-5037-T809070 | <input type="checkbox"/> uMic™ \$99.00 076-5027-T800070 |
| | <input type="checkbox"/> TVLink II \$99.00 076-0040-P511070 |

Step 5 Custom ear piece options

Shell style	xS 113/47	xP 127/57	xSP* 131/63	xSP†* 134/67
Hollow sleeve (hard) (HO) 062-0012-02	<input type="checkbox"/> L <input type="checkbox"/> R	NA	NA	NA
Solid sleeve (hard) (SO) 062-0012-02	<input type="checkbox"/> L <input type="checkbox"/> R	NA	NA	NA
Solid sleeve (soft) (S70) 062-0012-02	<input type="checkbox"/> L <input type="checkbox"/> R	NA	NA	NA
cShell (soft shell) (s60) 062-0014-02	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R	NA	NA
cShell (hard shell) 062-0014-02	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> L <input type="checkbox"/> R

*not available with Moxi Now

†only available on Moxi Dura

Receiver length

0 ☐ L ☐ R 1 ☐ L ☐ R 2 ☐ L ☐ R 3 ☐ L ☐ R

Domes

Open	<input type="checkbox"/> L <input type="checkbox"/> R	Small	<input type="checkbox"/> L <input type="checkbox"/> R	Medium	<input type="checkbox"/> L <input type="checkbox"/> R	Large
Closed	<input type="checkbox"/> L <input type="checkbox"/> R	Small	<input type="checkbox"/> L <input type="checkbox"/> R	Medium	<input type="checkbox"/> L <input type="checkbox"/> R	Large
Power	<input type="checkbox"/> L <input type="checkbox"/> R	Small	<input type="checkbox"/> L <input type="checkbox"/> R	Medium	<input type="checkbox"/> L <input type="checkbox"/> R	Large

cShell shell color/ Sleeve mold hard shell color (hard shell)

- ☐ Pink (26)
☐ Tan (14)
☐ Cocoa (22)
☐ Brown (28)
☐ Clear (21)*
☐ Blue/Red

cShell faceplate color (hard shell only)

- ☐ Pink (26)
☐ Tan (14)
☐ Cocoa (22)
☐ Brown (28)
☐ Clear (13)

Sleeve mold soft shell color

- ☐ Clear (21)*
☐ Trans. pink (T)
☐ Trans. brown (N)

Finish

☐ Gloss (HC) (standard) ☐ Satin (SA)

Canal length

<input type="checkbox"/> L <input type="checkbox"/> R Full length of impression	<input type="checkbox"/> L <input type="checkbox"/> R Long (RL)
<input type="checkbox"/> L <input type="checkbox"/> R Cut to mark	<input type="checkbox"/> L <input type="checkbox"/> R Medium (RM)
<input type="checkbox"/> L <input type="checkbox"/> R Extra long (RG)	<input type="checkbox"/> L <input type="checkbox"/> R Short (RS)

Venting

☐ L ☐ R Intellivent (audiogram required) (AO)

Select-A-Vent (S)

☐ L ☐ R Pressure vent 1.5 mm (S15)
☐ L ☐ R Small 1.8 mm (S18)
☐ L ☐ R Medium 2.5 mm (S25)
☐ L ☐ R Large 3.0 mm (S30)

☐ L ☐ R IROS vent (S, IA)
☐ L ☐ R As large as possible
☐ L ☐ R No vent (X)

Cavity vent (Y)*

☐ L ☐ R Pressure vent 1.5 mm (Y15)
☐ L ☐ R Small 1.8 mm (Y18)
☐ L ☐ R Medium 2.5 mm (Y25)
☐ L ☐ R Large 3.0 mm (Y30)

*Cavity vent default on hollow sleeve mold

Wax system

☐ L ☐ R None
☐ L ☐ R UH wax guard (CS)†
☐ L ☐ R Wax shield (ST)*

☐ L ☐ R Wax spring (WS)*
☐ L ☐ R Extended receiver tube (ER)*

*Only available on cShells

†Only available on sleeve molds

Shell options

☐ L ☐ R Removal filament (RF) ☐ L ☐ R Canal lock (CL) n/a with soft sleeve and cShell (soft). Same color as shell.

Step 6 Special instructions

Please send:

- | | |
|---|---------------------------------------|
| <input type="checkbox"/> Shipping labels | <input type="checkbox"/> Order forms |
| <input type="checkbox"/> Impression boxes | <input type="checkbox"/> Repair forms |

SPECIAL PROGRAMMING INSTRUCTIONS ON BACK

800.397.4442 | FAX 630.836.9770 | 444 Commerce Street | Aurora, Illinois 60504 | unitron.com/us

Special Programming Instructions

Set AAM to 100%	EXL04	<input type="checkbox"/>
NAL-NL1 fitting formula	FIT03	<input type="checkbox"/>
Add noise program	NOSP0	<input type="checkbox"/>
Increase wind reduction +1 Lev	WDB10	<input type="checkbox"/>
Decrease binaural MT gain - 12dB	BNP12	<input type="checkbox"/>
Set AAM to off	AAM00	<input type="checkbox"/>
Set AAM to 70%	AAM70	<input type="checkbox"/>
Increase gain in easy-t by 10dB	IET10	<input type="checkbox"/>
Deactivate program button	PBB00	<input type="checkbox"/>
Deactivate volume control	LVD00	<input type="checkbox"/>
Enable DuoLink for VC	LVD03	<input type="checkbox"/>
Short start up delay	SUD10	<input type="checkbox"/>
Long start up delay	SUD30	<input type="checkbox"/>
Turn data logging off	LOG00	<input type="checkbox"/>
Remove easy-t from program sequence	TCP00	<input type="checkbox"/>
Disable beep for program change	BPC00	<input type="checkbox"/>
Enable beep for warnings/alerts only	LBB01	<input type="checkbox"/>
Set push button to right louder left softer	PBB02	<input type="checkbox"/>
Disable all beep signals	LVD00	<input type="checkbox"/>
Turn frequency compression on		<input type="checkbox"/>
Add tinnitus program audio based		<input type="checkbox"/>
Add tinnitus program white noise		<input type="checkbox"/>
Add tinnitus program pink noise		<input type="checkbox"/>
Set push button to volume up and program change		<input type="checkbox"/>
Set push button to program change and volume down		<input type="checkbox"/>
Set push button to volume up		<input type="checkbox"/>
Set push button to volume down		<input type="checkbox"/>
Set volume control range to 20dB		<input type="checkbox"/>
CROS/BICROS		<input type="checkbox"/>

Moxi colors

