




Phonak Government Services RIC xReceiver and Coupling Guide

Fitting Range	xReceiver	Hearing Instrument		Coupling Options		
 Mild to moderately-severe hearing loss	Standard (xS) 46/111 (2cc) 	Phonak Audéo B  10 312 312T 13 R	Phonak Naida V  RIC	 Open Dome  Closed Dome  Power Dome	 SlimTip Hard  SlimTip Soft  cShell (Hard or Soft)	
 Mild to severe hearing loss	Power (xP) 57/124 (2cc) 	Phonak Audéo B  10 312 312T 13 R	Phonak Naida V  RIC	 Open Dome  Closed Dome  Power Dome	 SlimTip Hard  cShell (Hard or Soft)	
 Moderate to profound hearing loss	UltraPower (xUP) 66/130 (2cc) 	Phonak Audéo B  312 312T 13 R	Phonak Naida V  RIC	 cShell (Hard only)		

Phonak Government Services BTE Coupling Guide

Model	Earhook	SlimTube	Coupling Options
Phonak Bolero B  M P PR SP	 HE 10 (unfiltered) HE 10 680 (filtered)	 Soft SlimTube HE (0, 1, 2, 3)	 Open Dome Closed Dome Power Dome Phonak Free Field Phonak Skeleton Phonak Semi-Skeleton Phonak Canal Phonak Canal Lock SlimTip Hard SlimTip Soft Phonak Half Shell Phonak Carved Full Shell Phonak Standard Full Shell Phonak Helix Lock
Phonak Naida V  SP UP	 HE 11 (unfiltered) HE 11 680 (filtered)	 Power SlimTube II (00, 0, 1, 2, 3)	 Phonak Canal Phonak Skeleton Phonak Semi-Skeleton Phonak Helix Lock Phonak Half Shell Phonak Carved Full Shell Phonak Standard Full Shell Phonak Canal Lock
Phonak CROS B Phonak CROS II  312 13 312 13	 CROS Hook (13 model only)	 CROS SlimTube (0, 1, 2, 3)	 Open Dome Closed Dome Power Dome CROSTip